

Gamma SatPax[®] 670

COMPOSITION AND ATTRIBUTES

- Nonwoven cellulose/polyester
- Pre-wetted with EU Pharmacopoeia isopropyl alcohol (IPA) or Denatured Ethanol filtered to 0.22um and deionised water (DIH₂O)
- Z-Folded and wetted with measured volumes ensures consistently wetted wiper (60% saturation)
- Peel and Reseal label allows for single wipe use and prevents solvent from evaporation
- Flat packs allow for easy storage and reduce spill potential
- Non-standard mixtures are available
- Sterilisable by gamma irradiation
- Small pack sizes
- Double packed

APPLICATIONS

- Cleaning of exterior equipment in production environments
- Wiping and prepping of glassware and vials
- General wiping of surfaces

ADVANTAGES

- Convenient, ready to use format
- Absence of spray bottles reduces Volatile Organic Compound emissions
- Eliminates problems associated with filtration, sterilisation and mixing of isopropyl alcohol (IPA) or Denatured Ethanol.
- Optimises solution usage and reduces over-spray



VITAL STATISTICS (AVERAGE DATA IN UNSATURATED STATE)

Basis Weight (g/m ²)	68.9	Extractables:	
Caliper (µm)	247	H ₂ O (%)	0.090
Fibers (≥100µm/cm ²)	110	2-Propanol (%)	0.027
Particles (0.5µm) (biaxial shake)(x10 ³ /cm ²)	21		

ORDER INFORMATION

Product	IPA/DIH ₂ O	Item Number	Size /Metric	Shts/Pk	Pks/Cs
Gamma SatPax	70/30 %	GSPX670.001.0909.12	9x9"/23x23cm	25	12
Gamma SatPax	70/30 %	GSPX670.001.0909.40	9x9"/23x23cm	15	40

Product	DEB/DIH ₂ O	Item Number	Size /Metric	Shts/Pk	Pks/Cs
Gamma SatPax	70/30 %	GSPX670.011.0909.12	9x9"/23x23cm	25	12
Gamma SatPax	70/30 %	GSPX670.011.0909.40	9x9"/23x23cm	15	40