

SatPax[®] Lightweight Canister

This information is for the SatPax Lightweight Canister wiper pre-wetted with a solution of 70% pharmacopeia grade Isopropanol or Denatured Ethanol and 30% Deionised water. It has been developed for pharmaceutical and medical applications, offering a low cost, convenient and consistent wiper for general use. The canister is available as either standard, ready for irradiation or Gamma irradiated.

COMPOSITION AND ATTRIBUTES

- 23 gsm acrylic/pulp wiper
- 20cm x 20cm wiper size
- 200 sheets per canister
- Each standard canister has an outer bag.
- Each Gamma irradiated canister is Detex labeled and double bagged.
- Pre-Saturated with 70% pharmacopeia grade IPA or Denatured Ethanol and 30% Deionised water
- Convenient and reliable dispenser
- Perforated roll for easy tearing
- Can be Gamma Irradiated

APPLICATIONS

- ISO Class 4-7, B & C
- Non critical equipment and surface cleaning
- Pharmaceutical and medical
- Gamma Irradiated version suitable for sterile applications



VITAL STATISTICS

Berkshire fully tests all the critical parameters of the wiper as part of our ISO certified quality system. The data from these tests are not yet fully concluded, and will be made available as soon as complete.

ORDER INFORMATION

Product	Solution	Item Number	Size /Metric	Shts/Can	Can/Cs
SatPax Canister	70%IPA/30%DIW	SPXC23.01.10	8x8"/20x20cm	200	10
Ready for Irradiation	70%IPA/30%DIW	GSPXC23.01.10	8x8"/20x20cm	200	10
Gamma SatPax Canister	70%IPA/30%DIW	GSPXC23.01.10.ir	8x8"/20x20cm	200	10
SatPax Canister	70%DEB/30%DIW	SPXC23.11.10	8x8"/20x20cm	200	10
Ready for Irradiation	70%DEB/30%DIW	GSPXC23.11.10	8x8"/20x20cm	200	10
Gamma SatPax Canister	70%DEB/30%DIW	GSPXC23.11.10.ir	8x8"/20x20cm	200	10